

SHIELD^{me}™

DISPOSABLE ISOLATION GOWN

SDM-DGWN

CE



SHIELD^{me}™ NON-SURGICAL ISOLATION GOWN (WITH BELT) is intended to protect health care personnel and patients from the transfer of microorganisms, body fluids (Fluid repellent) and particulate material. Single-use and disposable. Latex Free.

Gown made of non-woven fabrics using spun bonded polypropylene spunbond non-woven, 40 gsm, 100 % new chip, coated with polyethylene film (28 gsm PP + 12 gsm PE), and latex free material with back closure, waste and neckties, featured sleeves closed with a belt on the right side. Liquid resistant due to polypropylene spun bond non-woven coated with Polyethylene film

CE-label according to EU Directive for Medical equipment 93/42/EEC, Annex VII, class 1, non-sterile.

SPECIFICATIONS:

Color	: Blue/Yellow
Gown Length	: 110 Cm Regular and 115 Cm Ex-Large
Storage	: Store under ambient temperature and dry condition, keep away from Direct Sunlight Durability unopened 5 years
Waste Management	: Polybag, combustion in incineration plants only creates carbon dioxide and water, which directly forms part of the natural circulation. Cardboard is combusted or recycled
Packing	: 10 pcs/Bag : 5 Bags/Carton (50 Gown / Carton)

**We reserve the right to change specifications of our product without prior notice.*

info@shieldmeglobal.com

www.shieldmeglobal.com