



# SHIELD<sup>me</sup>™

## FACE MASK

CE

**SHIELD<sup>me</sup>™ KN95 FILTERING FACE PIECE RESPIRATORS** (FFR), which are also called Disposable respirators/Face Mask, are subject to various regulatory standards around the world. KN95 respirators are commonly recommended by health authorities during pandemic or emergency situations. This Face Mask respirators certified and manufactured as per International Standard EN 149 : 2001 + A1 : 2009.

Manufactured environmental non-woven fabric, antibacterial meltblown filter material, 4-layer material, double-layer antibacterial meltblown cloth

### SPECIFICATIONS:

Filtering Efficiency:	96%
Aerosol Type:	Non-Oil
Color:	White
Size:	Universal Exhalation
Valve:	No Strap
Attachment:	Stapled
Adjustable Nosit Clips:	Yes
Standards:	KN95 (GB2626-2006)
Certifications:	FDA & CE Certified
Shelf Life:	3 years

## APPLICATIONS:

### **SHIELD<sup>me</sup><sup>TM</sup> KN95 FILTERING FACE PIECE RESPIRATORS/FACE MASK**

are commonly used for filtering non-oil-based particles such as those resulting from wildfires, PM 2.5 air pollution, volcanic eruptions, or bio-aerosols (e.g. viruses).



*\*We reserve the right to change specifications of our product without prior notice.*