

# SHIELDme™

## NON-COTACT IR THERMOMETER



CE

**SHIELDme™ NON-COTACT INFRARED THERMOMETER** is a medical device for body / forehead temperature checking, super fast one-second reading with high measurement accuracy.

### FEATURES:

With its innovative no-touch technology, this infrared thermometer is suitable for family members of all ages. It has smart built-in fever alarm system.

### SPECIFICATION PARAMETER:

Product Name	Non-Contact IR Thermometer
Product Model	SHIELDme
Body Temperature Range	34.0°C ~ 42.9 °C / 93 - 109 °F
Object Temperature Range	0°C - 100°C
Display Resolution	±0.2°C (0.2°F)
Measurement Distance	1 ~ 5 cm
Measurement Time	0.5 s
Automatic Power-off Function	15 s
Relative Humidity	≤85%
Storage Condition	-20°C ~ 50°C
Power Supply	DC 3.0V (2*AAA batteries)
Item Size	116 * 43 * 152mm (L * W * H)