

SHIELD^{me}TM

TOWER

SHIELD^{me}TM priority is your safety, **SHIELD^{me}TM** tower provides both high level disinfection and sanitization with the Automatic Hand Disinfecting Dispenser and body temperature detection through Touchless Temperature Measurement to ensure early detection and control of infection.

SHIELD^{me}TM tower can be used in all types of buildings, mainly where large people gather where it can be located in the entrance lobby, for everyone to sanitize upon entering the building, and at the same time making sure that whoever enters has a regular body temperature.

Height of the tower: 1790 mm
(without base height)

Width: 25 cm

Base dim: 400mm X 300mm X
100mm height

Depth of the tower: 4cm.

Materials: electrogalvanized
sheet

Thickness 1.2 mm.

Powder coated

Ral no. 9010

Off white



**Automatic
Hand Disinfecting
Dispenser**
(SDME AHDD1000)

- ✔ Kills 99.99% of viruses and bacteria
- ✔ Contains no alcohol
- ✔ Non flammable
- ✔ Non irritant
- ✔ Safe for humans
- ✔ Ecofriendly



**Touchless
Temperature
Measurement**
(SME IRTM03)

- ✔ Accurate
- ✔ Abnormal automatic alarm
- ✔ Automatic measurement
- ✔ Standby battery
- ✔ Operating temperature 10°C-40°C



**We reserve the right to change specifications of our product without prior notice.*

info@shieldmeglobal.com

www.shieldmeglobal.com